LCD HANDHELD DIGITAL MICROSCOPE

INSTRUCTION MANUAL

MODEL #44310

ENGLISH

Congratulations on your Celestron microscope purchase.

Your new microscope is a precision optical instrument made of high quality materials to ensure durability and long life, and is designed to give you a lifetime of pleasure with minimal maintenance.

Before attempting to use your microscope, please read through the instructions to familiarize yourself with the functions and operations to make usage easier for you. See the microscope diagram to locate the parts discussed in this manual.

The microscope provides optical powers of 1 to 3.7x and 54x on the LCD screen. It is ideally suited for examining solid objects such as coins, stamps, rocks, relics, insects, plants, skin, gems, circuit boards, various materials and many other objects. You can also examine some types of specimen slides at low and high powers.

You can observe the magnified images, capture video or take snapshots. Your images will be stored using the internal flash memory or with an optional SD card. You can upload your saved images to your computer with the supplied USB cable.

Note: This product was designed and intended for ages 13 and older!

DIAGRAM INFORMATION

1. 2.4" (60 mm) LCD Screen
2. "M" Button – Mode/Menu Settings
3. Power Button – On/Off
4. Snapshot/Video Button & OK/Confirm
5. Battery Cover
6. Up Button – Up and LEDs – 3 on/6 on/off
7. Down Button – Down & Digital Zoom
8. Change Power and Focus Ring
9. LED Cover
10. SD Card Slot
11. Ribs for Nonslip Handling
12. USB Port

SPECIFICATIONS

LCD Screen	2.4 inch /60 mm TFT 320 x 240 resolution
Sensor Resolution / Interpolation	3MP CMOS Sensor 2048x1536 / digital up to 12MP
Optical Power (Magnification)	1 to 3.7x and 54x with digital zoom up to 216x
Illumination	6 White LEDs
Digital Zoom	4x
Internal Storage	128MB Flash Memory (3MP – est. 250 Snapshots & 5min. Video)
SD Cards Supported	Up to 16GB
Field of View	26 mm @ 2x, 14 mm @ 3.7x, 2.5 mm @ 54x
USB Cable	2.0
Image Format	Snapshots: JPG — Video: AVI
Power Requirements	2AA Batteries (user supplied) 0.75W
Weight / Dimensions	6.6 oz (187 g) / 3.8" x 2.5" x 4.2" (96 x 63 x 107mm)
Operating Temperature	0°C to 40°C (32° F to 104° F)

GETTING STARTED

BATTERY INSTALLATION

To use your microscope install 2AA batteries which provides the power to operate it indoors or outdoors. Slide the battery cover on the bottom of the microscope (5) to open, and insert the batteries according to the battery polarity as shown. After the batteries are installed, close the battery cover. Battery condition will be shown in the top left of the LCD in viewing and video modes, and "Low Battery Condition" will be shown in red on the screen to notify you.

Note: When not using the microscope, you will want to turn it off to preserve battery life. Typical battery life is 1.5 hours using 3 LEDs and 0.75 hours using all 6 LEDs. For most applications, 3 LEDs are all that are needed for proper illumination.

Also, if the microscope will not be used for a few weeks or more, remove the batteries.

SET UP FUNCTIONS

- A. Viewing/Snapshot Mode** — “A” shown in the top right of the LCD screen. Press and hold for a couple of seconds the “M” button (2) to enter the Capture (Snapshot) and Set up functions. To change from “Capture” to “Set up” press and hold the “Down Button” (7) for a few seconds. To change from “Set up” to “Capture” press and hold the “Up Button” (6) for a few seconds. All settings are easy to do and very intuitive. In each function you press the “Up Button” or “Down Button” to navigate each choice, press the “Snapshot” button (4) to OK or confirm your choice.
- 1. Set Up Mode** — main set up items are the Language (choose among English, French, German, Italian, Spanish, Portuguese, Simplified Chinese, Traditional Chinese, Japanese) and Auto Off (time when the system turns itself off when not used for a specific time).
 - 2. Capture Mode** — main set up items are the Size (resolution settings from the sensor resolution of 3MP to lower and higher values through digital interpolation from 0.3MP up to 12MP), EV Value (Exposure Value) for increasing and decreasing the brightness level, Color choices and Time Stamp. You can exit back to the Viewing Mode by pressing the “Down Button.”
- B. Video Mode** — a Video Camera icon is shown in the top right of the LCD screen. From the Viewing/Snapshot Mode, press the “M” button to enter the Video Mode. Next, press and hold the “M” button to enter the Video Output and Set up functions. To change from “Video Output” to “Set up”, press and hold the “Down Button” for a few seconds and to change from “Set up” to “Video Output”, press and hold the “Up Button” for a few seconds. All settings are easy to do and very intuitive. For each function you press the “Up Button” or “Down Button” to navigate each choice, press the “Snap Shot” button to OK or confirm your choice.
- 1. Video Output** — set the “Size” (video resolution where VGA is the default) and EV (as described above).
 - 2. Set Up Mode** – same as described under Viewing/Snapshot mode for languages and auto shut off time.
- C. Playback Mode** — review the snapshots and videos you have taken in this mode. You can enter into this mode by pressing the “M” button and rotating from Viewing/Snapshot to Video to Replay.
- 1. Play** — you can Delete, Copy to Card or Protect images. Operation is similar to the other modes.
 - 2. Set Up Mode** — same as Viewing/Snapshot and Video set up modes Use the “Up Button” and “Down Button” to navigate through your images and videos in thumbnail sizes. To enlarge the images, press the “Snapshot/Video Button”. To delete images, use the Play mode above.

OPERATING THE MICROSCOPE

You can use your microscope to view objects, take snapshot images or videos to save and later upload to your computer. Whichever you choose to do, you first need to turn the power on by pressing (on some units press and hold for a few seconds) the

“Power” button (3). A red LED will light up under the button and shortly the screen will light up as three LEDs will illuminate the screen. To change the LEDs (the default setting is 3 LEDs on), press the “Up” button for 6 LEDs on. Press the button again, and the LEDs are turned off. To turn off the microscope, press (on some units press and hold) the “Power” button, and the screen will show “Shutting Power Off”. To change from Viewing Objects/Taking Snapshots to Taking Videos to Review Images, press the “M” button to rotate to the desired function.

A. Viewing Objects — To view objects, use the Viewing/Snapshot mode.

1. **Icons seen on the LCD Screen** — in the top right, “A” indicating viewing mode. In the top left, shows battery condition. The bottom right shows the resolution selected as well as the number of images available (to take) at the chosen resolution. Also at the bottom right, is “M” or “SD” (or similar designations) to note if internal memory or an optional SD Card is being used. On the left side, digital zoom power is displayed (if using this function).

2. **Focusing and changing Power** — use the Focus Ring (8) to change power and focus to see a sharp image. You can rest the LED Cover (9) directly on/over an object to view. Always start with the low power (3.7x) by rotating the focus ring to the far left and move the ring slightly to get a sharp focus. The most enjoyable views may be seen at low power which has a wider field of view and a brighter image. To achieve an even lower power setting (varying from 3.7x down to 1x), raise the microscope above the object. As the microscope is raised higher and power decreases, refocusing may be necessary. To change to high power (54x), rotate the Focus Ring all the way to the right and turn the ring slightly to obtain a sharp focus. You must use the high power directly on or slightly above an object to obtain a sharp focus. For low power, the three LEDs should be adequate, but you may need six LEDs when using high power to illuminate the object. Remember that using all six LEDs reduces the life of the batteries. To achieve higher powers than the standard low and high power, you can use the digital zoom feature by pressing the “Down Button” in the Viewing or Video modes. Pressing once will go to 2x digital zoom and pressing again will go to 4x. Pressing the button again reverses the zoom function. This feature gives you additional benefit by allowing you to digitally magnify an image to make it more enjoyable on some objects.

B. Taking Snapshot Images — To take snapshots, you must be in the Viewing/Snapshot mode. Press the “Snapshot/Video Button”, and the images will then be stored in internal memory or on an SD Card if using one.

C. Taking Videos — To take videos, you need to be in the Video Mode (the video icon will be in the top right of the LCD screen). Once the image is in focus, press the “Snapshot/Video Button” (4) and recording will begin. To stop the recording, press the same button again. The time of the video as it is being recorded, is shown at the bottom of the LCD.

D. Using an SD Card — You can also take images (snapshot or video) by using an SD (Secure Digital) Card. Your SD card is inserted in the SD Card Slot (10) on the

back of the microscope. When you insert the SD card properly, the SD card icon will be visible at the bottom right of the LCD screen. The bottom right of the screen will also show the memory available on the SD card. When imaging, the images are saved on the SD card and not in the internal storage memory. Do not pull and insert the SD Card quickly or errors may occur.

Note: Inserting or removing an SD card while the LCD is on, may cause the LCD to shut down and/or could damage the SD card.

E. Transferring Images to your Computer — To transfer images to a PC or MAC, you need to have a free USB port and have an imaging program to use for snapshots and/or videos.

Note: Do not disconnect the USB cable while transferring images or damage may occur.

You can transfer images from the internal storage memory to your PC by using the supplied Mini USB Cable. The small plug end of the cable plugs into the USB Port (12) on the back of the microscope, and the large plug end of the cable plugs into your computer. If the connections are proper, you will see (on your microscope screen) "MSDC" or similar data/icon. Your computer will automatically recognize the new hardware. Next, you will choose which program on your computer you want to transfer the images to.

If you used an SD card to store your images, you can transfer them to your PC using the same method mentioned above. You can also remove the SD card and insert it into your computer (if there is an SD card slot available) to transfer the images. Either way, your PC will ask you to choose which program you would like to transfer the images to.

CARE, MAINTENANCE AND WARRANTY

Your microscope is a precision optical instrument and should be treated with care at all times. Follow these care and maintenance suggestions and your microscope will need very little maintenance throughout its lifetime.

- Store the microscope in a dry and clean place inside the pouch included.
- Be very careful if using your microscope in direct sunlight to prevent damage to the microscope or your eyes.
- Never point the microscope towards the Sun, or the camera may become damaged and cease to operate.
- Clean the outside surfaces with a moist cloth.
- Blow off dust with a camel's hair brush or air blower from the optical surfaces.
- Never disassemble or clean internal optical surfaces. This should be done by qualified technicians at the factory or other authorized repair facilities.
- Do not disassemble or take apart the microscope or damage can occur.
- When handling glass specimen slides, use care as the edges can be sharp.
- If the screen freezes up for any reason, power off and start again or remove and replace the batteries.

WARRANTY:

YOUR MICROSCOPE HAS A TWO YEAR LIMITED WARRANTY.

PLEASE VISIT THE CELESTRON WEBSITE FOR DETAILED
INFORMATION ON ALL CELESTRON MICROSCOPES AT
WWW.CELESTRON.COM

FCC Statement

This device complies with Part 15 of FCC Rules. Operation is subject to the following two conditions:

1. This device may not cause harmful interference, and
2. This device must accept any interference received, including interference that may cause undesired operation.

2835 Columbia Street • Torrance, CA 90503 U.S.A.
Telephone: 310.328.9560 • Fax: 310.212.5835

©2011 Celestron
All rights reserved. • Printed in China • 03-11

Designed and intended for
those 13 years of age and older.

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